

# Multi Plankton Sampler



## Attention:

Please note that the Multi Plankton Sampler can only be operated upon request. Request operations with the Multi Plankton Sampler to AWI-Logistics in sufficient time prior to the cruise.  
Please ensure (prior to the previous shipyard) the winch AWI-004 is installed.

## Summary

The Multi Plankton Sampler MPS 92B is used to collect siliceous and calcareous proto-zooplankton (radiolarians and foraminifera) from different water depth levels to get information about the ecology of specific species and species groups used for geochemical and paleobiological applications.

<b>Manufacturer</b>	Hydro Bios
<b>Model</b>	MPS 92
<b>Serial No.</b>	n/a
<b>Type</b>	net
<b>REGISTRY-Link</b>	<a href="#">REGISTRY (8442)</a>



## Contacts

Name	Institution	Role
Barbara Niehoff	Alfred Wegener Institute, Helmholtz Centre for Polar and Marine Research	Device Contact

## Components

## Documentation

- [MultiNet\\_Midi Data Sheet \(Factsheet, 280 kB\)](#)
- [HYDRO-BIOS \(Web site\)](#)