

Multi Plankton Sampler



Attention:

Please note that the Multi Plankton Sampler can only be operated upon request. Request operations with the Multi Plankton Sampler to AWI-Logistics in sufficient time prior to the cruise.
Please ensure (prior to the previous shipyard) the winch AWI-004 is installed.

Summary

The Multi Plankton Sampler MPS 92B is used to collect siliceous and calcareous proto-zooplankton (radiolarians and foraminifera) from different water depth levels to get information about the ecology of specific species and species groups used for geochemical and paleobiological applications.

Manufacturer	Hydro Bios
Model	MPS 92
Serial No.	n/a
Type	net
REGISTRY-Link	REGISTRY (8442)



Contacts

Name	Institution	Role
Barbara Niehoff	Alfred Wegener Institute, Helmholtz Centre for Polar and Marine Research	Device Contact

Components

Documentation

- [MultiNet_Midi Data Sheet \(Factsheet, 280 kB\)](#)
- [HYDRO-BIOS \(Web site\)](#)